



Adhesive Surgical Dressing

Ideal for coping with all light to moderate exudate



FACT BOX

Product Code:
MD1155

Product Category:
Woundcare.

Product Approvals:
CE marked.

Size:
7cm x 8cm.

CAUTION:
SINGLE USE ONLY.
Danger of infection if reused.

The Medporex Adhesive Surgical Dressing is suitable for use on a wide variety of low to moderate exuding wounds, including surgical wounds, cuts and grazes. The dressing is made with a non-woven pad and material, allowing wounds to breathe to promote faster healing.

FEATURES




- Suitable for low to moderate exuding wounds.
- Sterile and individually packed.
- Non-adherent, non-woven pad.
- Breathable.
- 50 pieces per pack.

DIRECTIONS:

1. Peel back the film.
2. Place the pad over the wound area, peel back the remainder of the film carefully and keep the fingertips at the edge of the pad to avoid contamination.
3. Once the film is removed, gently press the pad on the wound area to secure placement.

STORAGE:

Store in a clean dry place with relative humidity.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid