



LIFESENSE®

A2 Plus Blood Pressure Monitor



FACT BOX

Product Code:
MD1803

Product Category:
Diagnostics

Product Approvals:
Medically certified under
EN ISO 13485:2016

Contents:
x1 Blood Pressure Monitor
x 1 Monitor Cuff (22~32cm)

Application Scope:
This device is to be used to measure the diastolic pressure, systolic pressure & pulse rate of adults by oscillometric method for diagnostic reference.



Scan above to download our free
MEDICARE LIFESENSE APP
Measure, Track and Share your Results



The Medicare LifeSense® A2 Plus Bluetooth Blood Pressure Monitor is easy to use, lightweight and portable, making it convenient to check your blood pressure from the comfort of your own home. It uses a comfortable upper arm cuff and presents your systolic, diastolic and pulse rate results on the easy-to-read digital display. The readings are also transferred by Bluetooth to the Medicare LifeSense® App on your smart device.

FEATURES

- Bluetooth enabled
- Clinically proven accuracy
- Irregular heartbeat detection
- WHO classification indicator
- Cuff circumference: 22~36cm
- Large text & soft touch buttons
- Accommodates 2 user profiles and holds 99 recordings per user
- Large easy-to-read display (60mm x 77mm)

WARNINGS/PRECAUTIONS:

Keep the device out of reach of children under 12 and people who cannot express their intention.

If children between the age of 12 and 18 use the device they must be accompanied by an adult.

Do not wrap the cuff over a wound.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid



LIFESENSE®

A2 Plus Blood Pressure Monitor




MEASURE, TRACK & SHARE YOUR RESULTS

The Medicare LifeSense® App saves all your blood pressure readings and creates a comprehensive chart to help create a picture of your long-term progress.

Compare your previous numbers or share your information with your healthcare professional via email from the app.

The app allows you to set reminders to check your blood pressure. You can also access all your information from the Medicare LifeSense® online portal at app.flemingmedical.ie.



1 Apply the cuff to your left arm with your palm facing up. Fasten the cuff around your arm, placed 2-3cm above the bend in your elbow.



2 While the device is turned off, press the "Start/Stop" button. The cuff around your arm will inflate automatically and stop when the measurement has been completed.



3 Once the measurement has been completed, your results will be displayed on the screen of the device and also synced to the app.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid