



Silicone Earplugs

Protects against water
and irritating sound



FACT BOX

Product Code:
MD71

Product Category:
Home Health

Contents:
3 pairs of Foam Silicone Earplugs

Danger:
The earplugs may be worn constantly in noisy environments.
Never use for diving.
Do not push the plug deep into the ear canal.
Keep away from heat and direct light.
Keep out of the reach of children.
Consult doctor before using if you have or have had ear problems (infections or perforated eardrum).
The product may be adversely affected by certain chemical substances.
Failure to comply with the above instructions for use may impair the efficiency of earplugs.
For further information contact Fleming Medical.

Medicare® Silicone Earplugs are designed to prevent disturbances in the hearing conduct caused by excessive noise or water. They are gentle, comfortable and adapt perfectly to the shape of the ear, offering excellent protection against external agents such as dust and wind.

FEATURES:

- One size fits all
- 3 pairs included
- Protects against water and irritating sound

DIRECTIONS:

1. Shape the ear protectors into a cone.
2. Place in front of the ear canal and gently press. Use entire piece per ear, do not split it.
3. Do not insert deep into the ear canal.

CONTRAINDICATIONS AND PRECAUTIONS

Keep out of the reach of children.
Do not push the plug deep into the ear canal.
Sudden removal may damage eardrum.
Never use for diving.
Keep away from direct light and heat.



@Medicarefirstaid



@Medicarefirstaid



@Medicare1staid